

Precautions

Benik's products DO NOT contain latex.

- Allergic reactions generally appear within the first few hours. Reactions that occur days or weeks later can often be attributed to heat, moisture or hygiene issues. Launder product regularly to avoid this issue. When cleaning product, do not deviate from recommended cleaning process.
- If unusual swelling, skin discoloration, rash or discomfort occur, discontinue use and schedule a follow up visit.
- Care should be taken to ensure that the support does not fit so tight that circulation is impaired.
- Not for people who have an allergic reaction to neoprene, or who have a tendency to develop dermatitis.
- Not for use over open wounds and infected areas.

Care and Cleaning

Wash frequently to remove body oils and other soil.

Hand wash in warm water and mild soap. Rinse thoroughly and air dry. Do not wring. Be sure support is completely dry before reapplying it.

Do not use ointments or oils under the support.

Do not dry clean or expose to heat or direct sunlight for extended periods of time.

Regular hygiene of the skin will prevent the inside of the product from becoming excessively soiled.

If the shape of the splint becomes deformed at any time, contact us for a follow up visit and we will reform the splint.

